

STERILIZATION POUCH AND SELF SEALING POUCHES

CHARACTERS AS FOLLOWS:

- 1- Pouch constructed from high grade medical Kraft paper and multilayer transparent film.
- 2- Strong seal strength.
- 3- Just like sterilization reel, fiber free opening.
- 4- Water based, nontoxic process indicators for steam and E.O gas are imprinted.
- 5- Excellent peel ability; quickly, Self Adhesive strip sealed.
- 6- Always keep the pouches clean and tidy.
- 7- Packed in boxes of 200 pcs.



SIZES OF OUR STERILIZATION FLAT POUCHES	
CODE	SIZES
IS 0520	50mm x 200mm
IS 7520	75mm x 200mm
IS 7525	75mm x 250mm
IS 7535	75mm x 350mm
IS 7540	75mm x 400mm
IS 1020	100mm x 200mm
IS 1025	100mm x 250mm
IS 1030	100mm x 300mm
IS 1035	100mm x 350mm
IS 1040	100mm x 400mm
IS 1045	100mm x 450mm
IS 1520	150mm x 200mm
IS 1525	150mm x 250mm
IS 1530	150mm x 300mm
IS 1535	150mm x 350mm
IS 1540	150mm x 400mm
IS 2025	200mm x 250mm
IS 2035	200mm x 350mm
IS 2040	200mm x 400mm
IS 2045	200mm x 450mm
IS 2535	250mm x 350mm
IS 2540	250mm x 400mm
IS 2545	250mm x 450mm
IS 2848	280mm x 480mm
IS 3040	300mm x 400mm
IS 3545	350mm x 450mm
IS 4065	400mm x 650mm

SIZES OF OUR STERILIZATION GUSSETED POUCHES	
CODE	SIZES



**IBN SINA MEDICAL
FACTOR**

Tel; 00973-17-731677, Fax; 00973-17-731675, web sit: WWW.IBNSINA-MEDICAL.COM

Email; Info@IBNSINA-MEDICAL.com, Sales@IBNSINA-MEDICAL.com

IS 7525	75mm x 250mm
IS 1020	100mm x 200mm
IS 1030	100mm x 300mm
IS 1045	100mm x 450mm
IS 1520	150mm x 200mm
IS 1530	150mm x 300mm
IS 1540	150mm x 400mm
IS 2039	200mm x 390mm
IS 2045	200mm x 450mm
IS 2050	200mm x 450mm
IS 2535	250mm x 350mm
IS 2545	250mm x 450mm
IS 3035	300mm x 350mm
IS 3045	300mm x 450mm
IS 4035	400mm x 350mm
IS 4045	400mm x 450mm