

STERILIZATION INDICATOR STRIP



Offers a complete line of color changing chemical indicator strips for the STEAM, EO, FORM, and PLASMA sterilization processes. Exposure to sterilizing conditions is confirmed conveniently through a distinct visual color change of the indicator ink on the strip. Long indicator strips come perforated allowing the use of either a full or half strip for differently sized packs. This chemical indicator monitors multiple sterilization parameters to assure sterilization conditions have been met..



Product Name	Color Change	QTY	mm x mm x mm MEAS	N.W. (Kg)	G.W. (Kg)	Operate Condition
STEAM	Yellow to Black	250pcs/box 8boxes/carton	152 x 90 x 170	1.8	1.4	121 C, 18min 132 C, 6min
EO	Red to Green	250pcs/box 8boxes/carton	152 x 90 x 170	1.8	1.4	(temperature) 54C (humidity)60% (time)90
FORM GAS	Red to Yellow	250pcs/box 8boxes/carton	152 x 90 x 170	1.8	1.4	(temperature) 75C (time) 1 h
PLASMA	Blue to Red	250pcs/box 8boxes/carton	152 x 90 x 170	1.8	1.4	30% (H ₂ O ₂), 50C, (Humidity) 95%



IBN SINA MEDICAL FACTORY

Tel; 00973-17-731677, Fax; 00973-17-731675, web sit: WWW.IBNSINA-MEDICAL.COM

Email; Info@IBNSINA-MEDICAL.com, Sales@IBNSINA-MEDICAL.com