

NEEDLES

NEEDLES' SIZES:

All sizes of Blister packed needles (16G ~ 33G), Disposable, sterile, blister, (100PCS\Box).



NEEDLES' PRODUCTION PROCESS



Needles are produced of steel which is heated till it's molten, then drawn through a die designed to meet the size requirements of the needle.

As it moves along the production line, the steel is further formed and rolled into a continuous, hollow wire.

The wire is appropriately cut to form the needle. Some needles are significantly more complex and are produced directly from a die casting. Other metal components of the needles are also produced in this manner.

PACKING:

Our products packed in blister bag consist of medical paper and film sealing together, printed with our logo, instructions, sizes, and sterilization date.



IBN SINA MEDICAL FACTORY

Tel; 00973-17-731677, Fax; 00973-17-731675

Email; Info@IBNSINA-MEDICAL.com, Sales@IBNSINA-MEDICAL.com